




CE
EN374
ISO13485

Single Use
Ambidextrous

Powdered


Disposable Vinyl Exam Glove 4.5g

Disposable Vinyl Exam Glove 4.5g

Powdered

Single Use
Ambidextrous

CE
EN374
ISO13485



Disposable Vinyl Exam Glove 4.5g

REF VM4540PD

L
100PCS

Disposable Vinyl Exam Glove 4.5g

REF VM4540PD

L
100PCS

Disposable Vinyl Exam Glove 4.5g

Powdered

Single Use
Ambidextrous

CE
EN374
ISO13485

Storage

Store in cool dry place.
Do not store gloves where temperatures may rise above 40 C(104 F).
Do not store in highly humid or damp areas.
Shield from direct sunlight, fluorescent lighting, x-rays, moisture and ozone.

Caution

Vinyl gloves are not intended for use as a chemical barrier.
Do not use these gloves for protection from hazardous chemicals.

LOT

ZH2208CS05



202208



202708

Manufacturer
Zhonghong Pulin Medical Products Co.,Ltd.
Address: West Industrial Park, Luannan County,
Tangshan City, Hebei Province, 063500 China
Phone: 86-0315-4167691

Made In China



AQL
1.5

