

## STAENDERMED - V

Velcro Supported Holders

Holds the canula with the agency of velcro (V) lock, 3 different models

### V1 Single piece mid-supported

The middle portion of the tube wraps, 2 sizes.

### V2 Single piece side-supported

The part close to the tube-to-end wraps 3 sizes.

### V3 Two pieces model

One size.


20



- Cotton surfaces touching the skin absorbs the sweat.
- Does not cause allergy and irritation in case of long term use.
- Cannula can be placed and removed without moving the patient's neck

- by the usage of single piece side-supported and two pieces models.
- Three different types of Endotracheal Cannula Holders: Velcro supported, silicon supported and clemp supported.
- Gamma sterilized.

Intensive Care Products



## STAENDERMED - SC

Silicone Supported Holder Bands

Holder band fixes the cannula by gripping it with its silicone surface (S). There are 3 different models.

### S1 Single piece, supported at middle part model.

The middle part wraps and grips the cannula. There are 2 sizes.

### S2 Single piece, supported by side part model.

The part close to one end wraps and grips the cannula. There are 3 sizes.

### S3 TP Two pieces model

There is only one standard single size.

### Staendermed Ordering Information

Model	Size (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code
V1	55	1	25	11 x 20 x 25	330 gr.	ETH 101
	65	1	25	11 x 20 x 25	345 gr.	ETH 102
	45	1	25	11 x 20 x 25	330 gr.	ETH 201
V2	55	1	25	11 x 20 x 25	345 gr.	ETH 202
	65	1	25	11 x 20 x 25	330 gr.	ETH 203
V3 TP	55	1	25	11 x 20 x 25	330 gr.	ETH TP
S1	55	1	25	11 x 20 x 25	330 gr.	ETH-SC 101
	65	1	25	11 x 20 x 25	345 gr.	ETH-SC 102
	55	1	25	11 x 20 x 25	330 gr.	ETH-SC 201
S2	65	1	25	11 x 20 x 25	345 gr.	ETH-SC 202
	55	1	25	11 x 20 x 25	330 gr.	ETH-SC 203
S3 TP	50 + 28	1	25	11 x 20 x 25	330 gr.	ETH-SC TP
CLP	55	1	25	11 x 20 x 25	330 gr.	ETH CLP 101
	65	1	25	11 x 20 x 25	345 gr.	ETH CLP 102
CLP TP	50 + 25	1	25	11 x 20 x 25	330 gr.	ETH CLP TP