



THREE-WAY STOPCOCK
sterile
TD-39

att.no.3 pt.2 subpt.2.3
edition III
date 2018-01-26
change of day: 2021-12-22

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.

ADDRESS: ul. Pod Borem 18, 41-808 Zabrze, Poland, phone +48 32 376 07 00-02, fax. +48 32 370 38 08

INDICATION:

- Device designed for managing the flow of infusion fluids during administration
- Allows to switch the infusion flow path
- For blood pressure monitoring
- For direct drug / fluid supply

FEATURES:

- Made of polycarbonate
- Transparent design
- Clear optical identifier of open/close position (marked flow directions) operating in 360 degrees
- Consistent, straight flow channels and an optimal internal diameter ensure correct pressure measurement
- Rotating Luer Lock
- All inlets protected with stoppers
- Resistant to lipids
- Latex and DEHP-free
- Single-use
- Non-pyrogenic, non-toxic
- EO sterilized
- Individual packaging: paper-foil



REFERENCE

CATALOGUE NUMBER	COLOUR	INDIVIDUAL PACKAGING	INTERMEDIATE PACKAGING	BULK PACKAGING
KT-002	blue	paper-foil	50 pcs	20 x 50 pcs