

Sterile Urethral Catheters

NANCHANG YILI MEDICAL INSTRUMENT CO., LTD

Office: 101# Building, Greenland International Exp Center, Jiulong Road 1177, Honggutan District, Nanchang, Jiangxi, China

Factory: Jinxian Medical Industrial Park, Nanchang, Jiangxi, China

Tel: +86-791-88252326 Fax: +86-791-88251326 Post Code: 330046

E-mail: sales@yilimed.com

www.yili-medical.com




■ Sterile Urethral Catheters

- Silicone Hematuria Catheter

1. Usage: prostatectomy, prostate cancer surgery, bladder cancer resection and urinary surgery operation.
2. 50 ml balloon for quick compression hemostasis.
3. Dufour tip and couvelaire tip for choose which providing a quick suction away blood clot and tissue after operation.
4. **3 way**, one for urine drainage, one for balloon inflation and deflation, and the third way for bladder irrigation and drug injection.



Couvelaire Tip



Dufour Tip

EasyThru®



Balloon Size: 50ML



Balloon Size: 30ML

Type	Code	Size (FR)	Length (MM)	Balloon (ML)	Color
	Normal Tube Wall				
Dufour Tip	Lot nr.120-122 YL-SFC3D-18	18	400- 410	50	Red
	YL-SFC3D-20	20	400- 410	50	Yellow
	YL-SFC3D-22	22	400- 410	50	Purple
	YL-SFC3D-24	24	400- 410	50	Blue
Couvelaire Tip	YL-SFC3C-18	18	400- 410	50	Red
	YL-SFC3C-20	20	400- 410	50	Yellow
	YL-SFC3C-22	22	400- 410	50	Purple
	YL-SFC3C-24	24	400- 410	50	Blue