

Combination of epT.I.P.S.[®] Pipette Tips with Eppendorf Research[®] 3 neo Pipettes

Eppendorf Research [®] 3 neo	0.1 – 10 µL, dark gray	0.1 – 10 µL M, medium gray	0.1 – 20 µL, medium gray	0.5 – 20 µL L, light gray	2 – 20 µL, yellow	2 – 100 µL, yellow	2 – 200 µL, yellow	20 – 300 µL, orange	50 – 1,000 µL, blue	50 – 1,250 µL, green	50 – 1,250 µL L, dark green	0.1 – 2.5 mL, red	0.1 – 5 mL, violet	0.2 – 5 mL L, violet	0.5 – 10 mL, turquoise	0.5 – 10 mL L, turquoise
epT.I.P.S. [®] BioBased	■		■	■			■	■	■	■	■	■	■	■		
epT.I.P.S. [®] Standard	■		■	■			■	■	■	■	■	■	■	■	■	■
epT.I.P.S. [®] BioBased Reloads			■		■	■	■	■	■	■	■					
ep Dualfilter T.I.P.S. [®]	■	■		■	■	■	■	■	■	■	■	■	■	■		■
ep Dualfilter T.I.P.S. [®] BioBased Reloads		■		■	■	■	■	■	■	■	■					
ep Dualfilter T.I.P.S. [®] SealMax [®]	■			■	■	■	■	■	■	■	■					
ep Dualfilter T.I.P.S. [®] SealMax [®] BioBased Reloads				■	■	■	■	■	■	■	■					
Variable volume, 1-channel pipettes																
0.1 – 2.5 µL	■	■	●	●												
0.5 – 10 µL	■	●	■	■	●											
1 – 20 µL	■	◆	◆	●	■											
1 – 20 µL	■				■	■	■	●								
5 – 100 µL	■				◆	■	■	●								
10 – 200 µL	■				◆	◆	■	●								
15 – 300 µL	■				◆	◆	◆	■								
50 – 1,000 µL	■								■	●	●					
0.1 – 2 mL	■											■				
0.25 – 5 mL	■												■	●		
0.5 – 10 mL	■														■	●

■ Test tip (other tips may require adjustment) ● Full volume range ◆ Limited volume range

Your local distributor: www.eppendorf.com/contact
 Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com

Eppendorf[®], the Eppendorf Brand Design, epT.I.P.S.[®], ep Dualfilter T.I.P.S.[®], SealMax[®] and Eppendorf Research[®] are registered trademarks of Eppendorf SE, Germany. All rights reserved, including graphics and images. Copyright © 2025 by Eppendorf SE. Order No: AA01 050 120 – EN1/0825/MCP



For more information go to
www.eppendorf.link/eptips

Combination of epT.I.P.S.[®] Pipette Tips with Eppendorf Research[®] plus Pipettes

Eppendorf Research [®] plus/ Eppendorf Research [®] plus Move It [®]	0.1 – 10 µL, dark gray	0.1 – 10 µL M, medium gray	0.1 – 20 µL, medium gray	0.5 – 20 µL L, light gray	0.1 – 20 µL, light pink	2 – 20 µL, yellow	2 – 100 µL, yellow	5 – 100 µL, light yellow	2 – 200 µL, yellow	20 – 300 µL, orange	50 – 1,000 µL, blue	50 – 1,250 µL, green	50 – 1,250 µL L, dark green	0.1 – 2.5 mL, red	0.1 – 5 mL, violet	0.2 – 5 mL L, violet	0.5 – 10 mL, turquoise	0.5 – 10 mL L, turquoise
epT.I.P.S. [®] BioBased	■		■	■					■	■	■	■	■	■	■	■		
epT.I.P.S. [®] Standard	■		■	■	■			■	■	■	■	■	■	■	■	■	■	■
epT.I.P.S. [®] BioBased Reloads			■			■	■		■	■	■	■	■					
ep Dualfilter T.I.P.S. [®]	■	■		■	■	■	■	■	■	■	■			■	■	■		■
ep Dualfilter T.I.P.S. [®] BioBased Reloads	■	■		■		■	■		■	■	■		■					
ep Dualfilter T.I.P.S. [®] SealMax [®]	■			■			■		■	■	■							
ep Dualfilter T.I.P.S. [®] SealMax [®] BioBased Reloads				■			■		■	■	■							
Fixed volume, 1-channel pipettes																		
10 µL	■	●	■	■	●													
20 µL	■		●	■														
10 µL	■					■	■		■	●								
20 µL	■					■	■		■	●								
25 µL/50 µL/100 µL	■						■		■	●								
200 µL	■								■	●								
200 µL/250 µL/500 µL/1,000 µL	■										■	●	●					
Variable volume, 1-, 8- and 12-channel pipettes																		
0.1 – 2.5 µL	■	■	●	●														
0.5 – 10 µL	■	●	■	■	●													
2 – 20 µL	■	◆	◆	●	■													
2 – 20 µL	■					■	■		■	●								
10 – 100 µL	■					◆	■		■	●								
20 – 200 µL	■					◆	◆		■	●								
30 – 300 µL	■					◆	◆		◆	■								
100 – 1,000 µL	■										■	●	●					
120 – 1,200 µL	■											●	■					
0.25 – 2.5 mL	■													■				
0.5 – 5 mL	■														■	●		
1 – 10 mL	■																■	●
Variable volume, 16- and 24-channel pipettes																		
1 – 20 µL	■				■													
5 – 100 µL	■							■										
Move It[®], variable volume, 4-, 6-, 8- and 12-channel pipettes with adjustable tip spacing																		
1 – 20 µL	■				■													
5 – 100 µL	■							■										
30 – 300 µL	■					◆	◆		◆	■								
120 – 1,200 µL	■											●	■					

■ Test tip (other tips may require adjustment) ● Full volume range ◆ Limited volume range

Your local distributor: www.eppendorf.com/contact
 Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com

Eppendorf[®], the Eppendorf Brand Design, epT.I.P.S.[®], ep Dualfilter T.I.P.S.[®], SealMax[®], Eppendorf Research[®], and Move It[®] are registered trademarks of Eppendorf SE, Germany. All rights reserved, including graphics and images. Copyright © 2025 by Eppendorf SE. Order No: AA01 050 120 – EN1/0825/MCP.



For more information go to
www.eppendorf.link/eptips

Combination of epT.I.P.S.[®] Pipette Tips with Eppendorf Reference[®] 2 Pipettes

Eppendorf Reference [®] 2	0.1 – 10 µL, dark gray	0.1 – 10 µL M, medium gray	0.1 – 20 µL, medium gray	0.5 – 20 µL L, light gray	2 – 20 µL, yellow	2 – 100 µL, yellow	2 – 200 µL, yellow	20 – 300 µL, orange	50 – 1,000 µL, blue	50 – 1,250 µL, green	50 – 1,250 µL L, dark green	0.1 – 2.5 mL, red	0.1 – 5 mL, violet	0.2 – 5 mL L, violet	0.5 – 10 mL, turquoise	0.5 – 10 mL L, turquoise
epT.I.P.S. [®] BioBased	■		■	■			■	■	■	■	■	■	■	■		
epT.I.P.S. [®] Standard	■		■	■			■	■	■	■	■	■	■	■	■	■
epT.I.P.S. [®] BioBased Reloads			■		■	■	■	■	■	■	■					
ep Dualfilter T.I.P.S. [®]	■	■		■	■	■	■	■	■			■	■	■		
ep Dualfilter T.I.P.S. [®] BioBased Reloads		■		■	■	■	■	■	■		■					
ep Dualfilter T.I.P.S. [®] SealMax [®]	■			■		■	■	■	■							
ep Dualfilter T.I.P.S. [®] SealMax [®] BioBased Reloads				■		■	■	■	■							
Fixed volume, 1-channel pipettes																
1 µL/2 µL	■	■	●	●												
5 µL/10 µL	■	●	■	■	●											
20 µL	■		●	■												
10 µL	■								●							
20 µL	■				■	■	■	■	●							
25 µL/50 µL/100 µL	■						■	■	●							
200 µL	■						■	■	●							
200 µL/250 µL/500 µL/1,000 µL	■									■	●	●				
2 mL/2.5 mL	■												■			
Variable volume, 1-, 8- and 12-channel pipettes																
0.1 – 2.5 µL	■	■	●	●												
0.5 – 10 µL	■	●	■	■	●											
2 – 20 µL	■	◆	◆	●	■											
2 – 20 µL	■					◆	◆	■	●							
10 – 100 µL	■					◆	■	■	●							
20 – 200 µL	■					■	■	■	●							
30 – 300 µL	■					◆	◆	◆	■							
100 – 1,000 µL	■									■	●	●				
0.25 – 2.5 mL	■											■				
0.5 – 5 mL	■												■	●		
1 – 10 mL	■														■	●

■ Test tip (other tips may require adjustment) ● Full volume range ◆ Limited volume range

Your local distributor: www.eppendorf.com/contact
 Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com

Eppendorf[®], the Eppendorf Brand Design, epT.I.P.S.[®], ep Dualfilter T.I.P.S.[®], SealMax[®] and Eppendorf Reference[®] are registered trademarks of Eppendorf SE, Germany. All rights reserved, including graphics and images. Copyright © 2025 by Eppendorf SE. Order No: AA01 050 120 – EN1/0825/BE



For more information go to www.eppendorf.link/eptips

Combination of epT.I.P.S.[®] Pipette Tips with Eppendorf Xplorer[®]/Xplorer[®] plus Pipettes

	0.1 – 10 µL, dark gray	0.1 – 10 µL M, medium gray	0.1 – 20 µL, medium gray	0.5 – 20 µL L, light gray	0.1 – 20 µL, light pink	2 – 20 µL, yellow	2 – 100 µL, yellow	5 – 100 µL, light yellow	2 – 200 µL, yellow	20 – 300 µL, orange	50 – 1,000 µL, blue	50 – 1,250 µL, green	50 – 1,250 µL L, dark green	0.1 – 2.5 mL, red	0.1 – 5 mL, violet	0.2 – 5 mL L, violet	0.5 – 10 mL, turquoise	0.5 – 10 mL L, turquoise
Eppendorf Xplorer [®] / Eppendorf Xplorer [®] plus/ Eppendorf Xplorer [®] plus Move It [®]																		
epT.I.P.S. [®] BioBased	■		■	■					■	■	■	■	■	■	■	■		
epT.I.P.S. [®] Standard	■		■	■	■			■	■	■	■	■	■	■	■	■	■	■
epT.I.P.S. [®] BioBased Reloads			■				■		■	■	■	■	■					
ep Dualfilter T.I.P.S. [®]	■	■		■	■		■	■	■	■	■	■	■	■	■	■		■
ep Dualfilter T.I.P.S. [®] BioBased Reloads		■		■			■		■	■	■	■	■					
ep Dualfilter T.I.P.S. [®] SealMax [®]	■			■			■		■	■	■							
ep Dualfilter T.I.P.S. [®] SealMax [®] BioBased Reloads				■			■		■	■	■							
Variable volume, 1-, 8- and 12-channel pipettes																		
0.5 – 10 µL	■	●	■	■	●													
1 – 20 µL	■	◆	◆	●	■													
5 – 100 µL	■						◆	■		●								
10 – 200 µL	■						◆	◆	■	●								
15 – 300 µL	■						◆	◆	■	■								
50 – 1,000 µL	■										■	●	●					
50 – 1,200 µL	■										◆	■	●					
0.1 – 2.5 mL	■													■				
0.2 – 5 mL															■	●		
0.5 – 10 mL																■	●	●
Variable volume, 16- and 24-channel pipettes																		
1 – 20 µL	■				■													
5 – 100 µL	■							■										
Move It[®], variable volume, 4-, 6-, 8- and 12-channel pipettes with adjustable tip spacing																		
1 – 20 µL	■				■													
5 – 100 µL	■							■										
15 – 300 µL	■						◆	◆	■	◆	■							
50 – 1,200 µL	■											■	●					

■ Test tip (other tips may require adjustment) ● Full volume range ◆ Limited volume range



Your local distributor: www.eppendorf.com/contact
 Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com

Eppendorf[®], the Eppendorf Brand Design, epT.I.P.S.[®], ep Dualfilter T.I.P.S.[®], Eppendorf Xplorer[®] and Move It[®] are registered trademarks of Eppendorf SE, Germany. All rights reserved, including graphics and images. Copyright © 2025 by Eppendorf SE. Order No: AA01 050 120 – EN1/0825

For more information go to
www.eppendorf.link/eptips