

iCover

Peripheral balloon expandable
PTFE covered stent system

Designed to
cover all needs



www.ivascular.global

iVascular[®]
therapies for living

Able to treat the most tortuous arteries

iCover

Designed to cover all needs


Enhanced visibility

Unique stent with radiopaque markers to facilitate the implantation and the post-dilation

- **ePTFE covered stent with high flexibility**


Due to:

- Open cell design and alternated links
- Nested peaks to avoid strut-to-strut contact



Best in class post-expansion capacity ($\geq 2\text{mm}$)

High ePTFE porosity to ensure covered stent
post-expansion capacity



MAXIMUM POST-EXPANSION DIAMETERS							
STENT DIMENSIONS		STENT DIAMETERS (mm)					
		5	6	7	8	9	10
STENT LENGHT (mm)	17	≤10mm					≤12mm
	27						
	37						
	57						

Excellent navigability and easy to implant

iCover
Designed to cover all needs

Outstanding flexibility

of the crimped stent, avoiding the kinking during navigation

Small profiles:

6F introducer compatibility up to 8mm*

- **Crimped on an Oceanus 35 PTA balloon, offering:**

- Short deflation times
- Good pushability and trackability

INTRODUCER COMPATIBILITIES							
STENT DIMENSIONS		STENT DIAMETERS (mm)					
		5	6	7	8	9	10
STENT LENGHT (mm)	17	6F			7F		
	27						
	37						
	57						


* 6F up to 8x17 mm and 7F from 8x27mm up to 10 mm.

iVascular proprietary encapsulation technology

iCover
Designed to cover all needs



Specific ePTFE laser cutting technology




Encapsulation technology

iCover is a stent completely encapsulated with 2 ePTFE layers



Proprietary technology to encapsulate the stent into an inner and outer ePTFE layer

Ensures complete encapsulation avoiding delamination



• **Perfect sealing**

The laser follows the stent outline while cutting the contour

iCover features

- > Over the wire catheter (OTW)
- > Guidewire compatibility: **0.035"**
- > Catheter length: **80** and **140** cm
- > NP: 9 atm
- > RBP:
 - Ø 5-7mm: 15 atm
 - Ø 8mm: 14 atm
 - Ø 9-10mm: 13 atm
- > Maximum post-expansion:
 - Ø 5-8mm: ≤10mm
 - Ø 9-10mm: ≤ 12mm
- > CoverTech encapsulation technology
- > Introducer compatibility
6F (up to Ø8mm, L 17mm), **7F** (≥8mmØ, L27mm)
- > **3 radiopaque markers** at each stent end
- > Stent material: CoCr L605
- > Graft material: ePTFE
- > Recoil: ≤ 10 %
- > Shortening:
 - Length 17mm: ≤ 12 %
 - Length 57mm: ≤ 8 %

Product with CE mark, certified by Notified Body 0318

Effective catheter length (cm)	Stent diameter (mm)	Stent length (mm)			
		17	27	37	57
80	5	SPCBCC35080050017	SPCBCC35080050027	SPCBCC35080050037	SPCBCC35080050057
	6	SPCBCC35080060017	SPCBCC35080060027	SPCBCC35080060037	SPCBCC35080060057
	7	SPCBCC35080070017	SPCBCC35080070027	SPCBCC35080070037	SPCBCC35080070057
	8	SPCBCC35080080017	SPCBCC35080080027	SPCBCC35080080037	SPCBCC35080080057
	9		SPCBCC35080090027	SPCBCC35080090037	SPCBCC35080090057
	10		SPCBCC35080100027	SPCBCC35080100037	SPCBCC35080100057
140	5	SPCBCC35140050017	SPCBCC35140050027	SPCBCC35140050037	SPCBCC35140050057
	6	SPCBCC35140060017	SPCBCC35140060027	SPCBCC35140060037	SPCBCC35140060057
	7	SPCBCC35140070017	SPCBCC35140070027	SPCBCC35140070037	SPCBCC35140070057
	8	SPCBCC35140080017	SPCBCC35140080027	SPCBCC35140080037	SPCBCC35140080057
	9		SPCBCC35140090027	SPCBCC35140090037	SPCBCC35140090057
	10		SPCBCC35140100027	SPCBCC35140100037	SPCBCC35140100057

● Introducer compatibility: **6F** ● Introducer compatibility: **7F**

Distributed by:



Manufactured by:
Life Vascular Devices Biotech, S.L.
www.ivascular.global
info@ivascular.global