


Arterial



Venous

Arterial Tubing (red)						
No.	ERP	ERP	Old ERP	Descriptions	Qty	Mark
1	4060100001		04.010207	Cap of dialyser connector	1	
2	4070100008		04.021605	Dialyser connector (red)	1	
3	3030200008	4010100001	04.030107	Drip Chamber(OD:19-30mm, L:130mm)	1	
		4020200001	04.040307	Cover for Drip Chamber		
4	4040100003		04.062605	Clamp (red small)	1	
5	3032000008	4070600003	04.100205	Female luer lock Ø4.1	1	
		4060200001	04.112807	Cap of Female luerlock		
6	4070200001		04.080107	Pump Segment Connector Ø4.1	1	
7	4070200003		04.080205	Pump Segment Connector Ø4.1	1	
8	3032300007	4040100003	04.062605	Clamp (red small)	2	A
		4070600003	04.100205	Female luer lock Ø4.1		
		4060200001	04.112807	Cap of Female luerlock		
		3020600002		Branch Tubing: 2.2*4.1*100mm		
9	3032400007	4040100003	04.062605	Clamp (red small)	1	B
		4070600027	04.103105	Female luer lock Ø2.5		
		4060200001	04.112807	Cap of Female luerlock		
		3020600018		Heparin Line: 1.0*2.5*600mm		
		3020600001		Branch Tubing: 2.2*4.1*50mm		C
10	3030400001		04.140305	Access Port(red)	1	
11	4040200003		04.064305	Clamp (red big)	1	
12	3030500001		04.182705	Self-ejecting Patient Connector (red)	1	
13	3030900001	4070600003	04.100205	Female luer lock Ø4.1	1	
		4040100003	04.062605	Clamp (red small)		
		3030800001	04.221807	Suspended Spike		
		3020300010		Branch Tubing: 2.5*4.1*100mm		
14	3020400033			Main Tubing: 4.5*6.8*1000mm	2	
15	3020600008			Branch Tubing: 2.2*4.1*700mm	1	
16	3020400023			Main Tubing: 4.5*6.8*600mm	1	
17	3021100007			Pump Segment: 6.3*9.8*345mm	1	
18	3020400028			Main Tubing: 4.5*6.8*800mm	1	

Venous Tubing (blue)						
No.	ERP	ERP	Old ERP	Descriptions	Qty	Mark
1	4060100001		04.010207	Cap of dialyser connector	1	
2	4070100009		04.021704	Dialyser connector (blue)	1	
3	3032300042	4060200007	04.113607	Cap of Female luerlock	1	B
		3030300001	04.070907	Transducer Protector		
		4070600005	04.100304	Femalue luer lock Ø4.1		
		4040100005	04.062704	Clamp (blue small)		
		3020600005		Branch Tubing: 2.2*4.1*400mm		
4	3030200042	4020200001	04.040307	Cover for Drip Chamber	1	
		4010100001	04.030107	Drip Chamber(OD:19-30mm, L:130mm)		
		4030100001	04.050506	Filter		
5	3030400003		04.140404	Access Port (blue)	1	
6	3032300009	4060200001	04.112807	Cap of Female luerlock	1	A
		4070600005	04.100304	Femalue luer lock Ø4.1		
		4040100005	04.062704	Clamp (blue small)		
		3020600002		Branch Tubing: 2.2*4.1*100mm		
7	4040200005		04.064404	Clamp (blue big)	1	
8	3030500002		04.182804	Patient connector (blue)	1	
9	4071200002		04.230507	Recirculating connector	1	
10	3020400016			Main Tubing: 4.5*6.8*400mm	1	
11	3020400020			Main Tubing: 4.5*6.8*500mm	1	
12	3020400051			Main Tubing: 4.5*6.8*1800mm	1	

10	Rev.	Version:A.7Q	2021/10/30	lai	Prepared by :	Reviewed by :	Description:		
9	Rev.	Version:A.6Q	2020/11/10	liu	Name		<b>Tubing Sets for Hemodialysis</b>  <b>Internal ID:</b>  <b>BLU017E</b>  <b>Blood pathway</b>		
8	Rev.	Version:A.5Q	2020/1/10	lai	Date				
7	Rev.	Version:A.4Q	2019/6/14	hou	Approved by QA:				
6	Rev.	Version:A.3Q	2019/5/31	hou	Name				
5	Rev.	Version:A.2Q	2019/4/25	hou	Date		<b>Fresenius 4008</b>		
4	Rev.	Version:A.1 Q	2018/8/27	duan					
01	New Release		2016/9/30	hou			<b>STERILIZATION</b> <b>EO</b> <b>DEHP</b>		
Rev	Remarks	Version No.	Date	Name					
<b>CARTON SIZE</b>		<b>564*381*249mm</b>		<b>QTY/CTN</b>	<b>24sets</b>				

The customer information (客户信息)					
No. (序号)	Countries&Customer (国家客户)	Customer product code (客户产品代码)	Customer drawings model (客户图模型号)	Machine type (机型)	Note (备注)