


# medbar®

## Limb Holder




REF 238



CE 2292

STERILE EO



### Intended Use

Limb Holders are used for fastening the patient's hands and feet having a surgery or confined to bed who might hurt himself or environment unconsciously. Limb Holders are sterile disposable devices with single use.

### Indications

Limb Holders are used for fastening the patient's hands and feet where lying, to prevent a possible damage may be caused by patient unconsciously.

### Quality and Legistion

- CE marked according to 93/42/EEC Medical Device Directive and classified as Class I sterile medical device.
- Manufactured under EN ISO 13485:2012 Quality Management System Standard.
- Notified Body: UDEM Ltd. AŞ. (2292)
- EC Certification No: M.2016.106.7000

### Biocompatibility

- Biological evaluation of Limb Holders done according to ISO 10993-1.
- The products do not incorporate a substance of a human blood derivative, animal orginated tissues,phthalates, medicinal product, latex, radioactive materials and electromagnetic waves.

### Sterilization

- Limb Holders are sterilized with Ethylene oxide. Validated for ISO 11135 standard.



| REF NO | PRODUCT NAME | FEATURES      | STERILIZATION | PACKAGING UNITS  |
|--------|--------------|---------------|---------------|--|
| 23801  | Limb Holder  | • Adult Hand. | ETO           | 25 pcs / inner bag<br>250 pcs / outer box<br>Box Dimensions (cm): 43*37*56<br>Box Weight: 7,00 kg  |
| 23802  | Limb Holder  | • Adult Foot. | ETO           | 25 pcs / inner bag<br>250 pcs / outer box<br>Box Dimensions (cm): 43*37*56<br>Box Weight: 7,60 kg  |
| 23803  | Limb Holder  | • Kids.       | ETO           | 100 pcs / inner bag<br>500 pcs / outer box<br>Box Dimensions (cm): 39*43*41<br>Box Weight: 6,50 kg |
| 23804  | Limb Holder  | • Baby.       | ETO           | 100 pcs / inner bag<br>500 pcs / outer box<br>Box Dimensions (cm): 39*43*41<br>Box Weight: 6,50 kg |

