



## LB6810

### Closed Suction Catheter

- Sputum suction pipe sealing film produced by imported resin, strong toughness, feel is exquisite and conducive to health care workers
- Sputum suction tube body design scale mark line, facilitate throughout the sputum suction effect
- Closed suction phlegm tube adopts full closed type design, can be used to sync with breathing machine, at the time of sputum suction oxygen to sustainable, convenient and safe operation, to avoid cross-infection
- Specification: Fr6-Fr16
- 24 hours /72 hours

Cat.No.	Size (Fr)	Specification
LB681006	6	Tube diameter 2mm, length 450mm,pediatric
LB681008	8	Tube diameter 2.6mm, length 450mm,pediatric
LB681010	10	Tube diameter 3.3mm, length 600mm,Adult
LB681012	12	Tube diameter 4mm, length 600mm,Adult
LB681014	14	Tube diameter 4.7mm, length 600mm,Adult
LB681016	16	Tube diameter 5.3mm, length 600mm,Adult



## LB6811

### Closed Suction Catheter

- Sputum suction pipe sealing film produced by imported resin, strong toughness, feel is exquisite and conducive to health care workers
- Sputum suction tube body design scale mark line, facilitate throughout the sputum suction effect
- Closed suction phlegm tube adopts full closed type design, can be used to sync with breathing machine, at the time of sputum suction oxygen to sustainable, convenient and safe operation, to avoid cross-infection
- Specification: Fr6-Fr16
- 24 hours /72 hours

Cat.No.	Size (Fr)	Specification
LB681106	6	Tube diameter 2mm, length 450mm,pediatric
LB681108	8	Tube diameter 2.6mm, length 450mm,pediatric
LB681110	10	Tube diameter 3.3mm, length 600mm,Adult
LB681112	12	Tube diameter 4mm, length 600mm,Adult
LB681114	14	Tube diameter 4.7mm, length 600mm,Adult
LB681116	16	Tube diameter 5.3mm, length 600mm,Adult