



elastopor EYE
eye dressing, non-woven, with absorbent pad,
self-adhesive, sterile

att.no.3 pt.2 subpt.2.3

revision 2

date: 2016-06-03

TD-16

amendment of: 2023-02-01

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.

ADDRESS: ul. Pod Borem 18, 41-808 Zabrze, Poland, phone +48 32 376 07 00-02,06 fax. +48 32 370 38 08

INDICATION:

- covering an eye and its area after ophthalmic procedures and surgeries
- sight correction in strabismus therapy

FEATURES:

- self-adhesive eye dressing with absorbent made of hydrophobic non-woven material
- covered with acrylic glue
- micropores in the non-woven material structure provide an adequate rate of vapour permeability so correct gas exchange between the dressing and skin is guaranteed
- centrally situated absorbent pad absorbs moderate amounts of exudates and protects against external influences
- absorbent pad with polyethylene mesh does not stick to the wound
- elastic - comfortable for the patient, adjusts to the body shape and guarantees freedom of movement
- rounded edges prevent accidental removal of the dressing
- silicone paper backing facilitates precise, painless and efficient application
- available in various sizes and colours of non-woven and absorbent pad
- hypoallergenic
- sterile
- shelf life: 5 years



SIZE AND REFERENCE

CATALOGUE NUMBER	DIMENSIONS	SHAPE	COLOUR	INDIVIDUAL PACKAGING PAPER-PAPER	SALES PACKAGING BOX	TRANSPORT PACKAGING
803010	5,8 cm x 8,3 cm		absorbent pad - white	1 pc.	50 pcs	24 x 50 pcs
803011	6,5 cm x 9,5 cm		non-woven material - white	1 pc.	50 pcs	24 x 50 pcs
803212	5 cm x 7,5 cm		absorbent pad - white non-woven material - skin colour	1 pc.	50 pcs	24 x 50 pcs



Adhesive non-woven tape SOFTplast

TD-06-I.2.c-5.6b

edition I

date: 2021-05-26

amendment of: 2021-12-22

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.**ADDRESS:** ul. Pod Borem 18, 41-808 Zabrze, Poland; T: +48 32 376 07 00-02,06 F: +48 32 370 38 08; www.zarys.com**CLASSIFICATION**

class I medical device compliant with Medical Devices Regulation (EU) 2017/745 (MDR)

INTENDED PURPOSE

- device intended for stabilisation, protection and fixation of drains, catheters, medical tubes not requiring specialised fixation at their place of application as well as primary dressings (contact), non-adhesive, applied directly to wounds and skin damages

FEATURES

- made of non-woven fabric
- white colour
- wound on a roll
- covered with an acrylic adhesive applied by the continuant method
- can be divided lengthwise and crosswise without the use of scissors
- straight edge finishing
- good air and water vapour permeability
- good adhesion
- thin, flexible, conforms to the body contours, guarantees freedom of movement
- hypoallergenic
- single-use
- non-sterile
- shelf life: 3 years



CATALOGUE NUMBER	DIMENSIONS	INTERMEDIATE PACKAGING (SALES) BOX	BULK PACKAGING CARTON
SNW12505W	1,25 cm x 5 m	24 x 1 pc.	30 x 24 pcs
SNW12509W	1,25 cm x 9,14 m	24 x 1 pc.	30 x 24 pcs
SNW25005W	2,5 cm x 5 m	12 x 1 pc.	30 x 12 pcs
SNW25009W	2,5 cm x 9,14 m	12 x 1 pc.	30 x 12 pcs
SNW50005W	5 cm x 5 m	6 x 1 pc.	30 x 6 pcs
SNW50009W	5 cm x 9,14 m	6 x 1 pc.	30 x 6 pcs



elastopor IV
IV cannula dressing, non-woven,
self-adhesive, sterile

att.no.3 pt.2 subpt.2.3

revision 3

date: 2020-01-10

TD-09

amendment of: 2023-02-01

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.


ADDRESS: ul. Pod Borem 18, 41-808 Zabrze, Poland, phone +48 32 376 07 00-02,06 fax. +48 32 370 38 08

INDICATION:

- stabilization, protection and fixation of intravenous cannulas with a vertical port

FEATURES:

- self-adhesive dressing having a notch for the vertical port of the intravenous cannula
- made of hydrophobic non-woven material
- micropores in the non-woven material structure provide an adequate rate of moisture vapour permeability so correct gas exchange between the dressing and skin is guaranteed
- centrally situated absorbent pad protects the non-woven material from sticking to the injection site
- non-woven lining should be placed under cannula's wings, which prevents it from causing damage to the skin
- hypoallergenic
- elastic - adjusts to the body shape and guarantees freedom of movement
- easy to attach to the skin
- rounded edges prevent accidental removal of the dressing
- coated with a hypoallergenic acrylic adhesive without zinc oxide, applied to the nonwoven fabric by a continuous method
- silicone paper backing facilitates precise, painless and efficient application
- graphical instructions for use on individual and intermediate packaging
- EO sterilized
- shelf life: 5 years



SIZE AND REFERENCE

CATALOGUE NUMBER	DIMENSIONS	INDIVIDUAL PACKAGING PAPER-PAPER	SALES PACKAGING BOX	TRANSPORT PACKAGING
811013	6 cm x 8 cm	1 pc.	100 pcs	20 x 100 pcs
811012	5,1 cm x 7,6 cm	1 pc.	100 pcs	20 x 100 pcs
811023	6 cm x 8 cm	1 pc.	50 pcs	40 x 50 pcs
811022	5,1 cm x 7,6 cm	1 pc.	50 pcs	40 x 50 pcs



elastoDERM F-IV
IV cannula dressing, foil, with frame,
with U-incision, self-adhesive, sterile

TD-18

att.no.3 pt.2 subpt.2.3

revision 2

date: 2016-06-03

amendment of: 2023-02-01

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.

ADDRESS: ul. Pod Borem 18, 41-808 Zabrze, Poland, phone +48 32 376 07 00-02,06 fax. +48 32 370 38 08

INDICATION:

- fixing and protecting intravenous cannulas inserted into peripheral and central veins, peripheral entry into central veins
- fixation of pulmonary artery catheters, haemodialysis catheters and catheters used during anaesthetic and cardiac procedures

FEATURES:

- made of polyurethane film
- covered with hypoallergenic acrylic glue
- transparent – allows observation of the dressed place without removing the dressing
- vapour permeability guarantees correct gas exchange between the foil and skin and prevents it from excessive accumulation of moisture and maceration
- barrier against microorganisms and bacteria minimizing infection risk
- skin contact up to 7 days
- very thin and elastic – adjusts to the body shape and guarantees freedom of movement
- the application system – „frame” type – facilitates precise and painless application of the dressing
- waterproof
- adhesive
- hypoallergenic
- shelf life: 3 years



SIZE AND REFERENCE

CATALOGUE NUMBER	DIMENSIONS	INDIVIDUAL PACKAGING PAPER-PAPER	SALES PACKAGING BOX	TRANSPORT PACKAGING
812006	6 cm x 7 cm	1 pc.	100 pcs	30 x 100 pcs



elastopor STERIL
non-woven dressing with absorbent pad,
self-adhesive, sterile

TD-09

att.no.3 pt.2 subpt.2.3

revision 2

date: 2021-10-29

amendment of: 2023-02-01

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.

ADDRESS: ul. Pod Borem 18, 41-808 Zabrze, Poland, phone +48 32 376 07 00-02,06 fax. +48 32 370 38 08

INDICATION:

- treating wounds with a small and medium amount of exudate
- especially recommended for covering surgical wounds characterized by straight cuts and the application of fewer or more skin sutures in a straight line



FEATURES:

- self-adhesive dressing with absorbent pad
- made of hydrophobic non-woven material
- micropores in the non-woven material structure provide an adequate rate of moisture vapour permeability so correct gas exchange between the dressing and skin is guaranteed
- centrally placed absorbent pad - absorbs medium wound exudate
- absorbent pad coated with polyethylene mesh protects against external factors' influence and does not stick to the wound as well as facilitates removal of the dressing
- hypoallergenic
- elastic - comfortable for the patient, adjusts to the body shape and guarantees freedom of movement
- easy to attach to the skin
- rounded edges prevent accidental removal of the dressing
- covered with hypoallergenic acrylic glue without the content of zinc oxide, rubber, latex, evenly applied over the entire surface, with very good moisture resistance ensuring good stabilization of the dressing at the application site
- silicone paper backing facilitates precise, painless and efficient application
- EO sterilized
- shelf life: 5 years

SIZE AND REFERENCE

CATALOGUE NUMBER	DIMENSIONS	INDIVIDUAL PACKAGING PAPER-PAPER	SALES PACKAGING BOX	TRANSPORT PACKAGING
801001E	5 cm x 7,2 cm	1 pc.	100 pcs	24 x 100 pcs
801012E	6 cm x 10 cm	1 pc.	50 pcs	45 x 50 pcs
801002E	8 cm x 10 cm	1 pc.	30 pcs	36 x 30 pcs
801013E	10 cm x 10 cm	1 pc.	30 pcs	24 x 30 pcs
801021E	10 cm x 12 cm	1 pc.	30 pcs	30 x 30 pcs
801014E	10 cm x 15 cm	1 pc.	30 pcs	24 x 30 pcs
801003E	6 cm x 15 cm	1 pc.	50 pcs	24 x 50 pcs
801004E	8 cm x 15 cm	1 pc.	30 pcs	24 x 30 pcs
801005E	8 cm x 20 cm	1 pc.	30 pcs	18 x 30 pcs
801015E	10 cm x 20 cm	1 pc.	25 pcs	18 x 25 pcs
801016E	10 cm x 25 cm	1 pc.	25 pcs	18 x 25 pcs
801017E	10 cm x 30 cm	1 pc.	25 pcs	12 x 25 pcs
801018E	10 cm x 35 cm	1 pc.	25 pcs	12 x 25 pcs

Marketing material intended for the health care professionals.

Prepared by: W.W. – Product Manager ZARYS International Group sp. z o.o. sp.k.

No part of this document may be changed or copied without the permission of ZARYS International Group sp.z o.o. sp.k.

page 1 / 1

rev.2023.01



elastoLUMENAL S
eye pad, superabsorbent, multi-layer, sterile

TD-16

att.no.3 pt.2 subpt.2.3

revision 2

date: 2016-06-03

amendment of: 2023-02-01

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.

ADDRESS: ul. Pod Borem 18, 41-808 Zabrze, Poland, phone +48 32 376 07 00-02,06 fax. +48 32 370 38 08

INDICATION:

- covering an eye and its area after ophthalmic procedures and surgeries
- sight correction in strabismus therapy

FEATURES:

- composed of external non-woven layers and internal multi-layer cellulose pulp
- needs to be fixed with adhesive tape or in different way, e.g. with tubular net dressing
- hypoallergenic
- sterile
- shelf life: 5 years



SIZE AND REFERENCE

CATALOGUE NUMBER	DIMENSIONS	INDIVIDUAL PACKAGING PAPER-PAPER	SALES PACKAGING BOX	TRANSPORT PACKAGING
803310	5,5 cm x 7,5 cm	1 pc.	50 pcs	24 x 50 pcs