

# Specifications

Resolution*	854×480 pixel (NAVI-60 )
Image Light	10-25 Lumens 4 level adjustability
Color	5 colors
Image Size	3 Sizes
Vein depth detection	Yes
Screen*	Yes
Charge during operation	Yes
Size	224 (W) ×68(H)×64(D)mm
Weight	500g

\*NAVI-60 only



NAVI-FS Fixed Support



NAVI-MS Mobile Support



VeinNavi®  
Vein Illuminator  
NAVI-60/30



## MEDCAPTAIN MEDICAL TECHNOLOGY CO., LTD.




Add: 12th Floor, Baiwang Research Building, No.5158 Shahe West Road, Xili, Nanshan, 518055 Shenzhen, Guangdong, PEOPLE'S REPUBLIC OF CHINA  
Tel: +86 755 26953369 E-mail: info@medcaptain.com  
Fax: +86 755 26001651 Website: www.medcaptain.com  
MEDCAPTAIN is a trade mark of Medcaptain Medical Technology Co., Ltd.  
Specifications are subject to change without prior notice.  
©2019 Medcaptain Medical Technology Co., Ltd. All rights reserved.  
EN-NAVI-60/30-2P-Version 2.1

DISTRIBUTOR:

**NAVI-60 Vein Illuminator** provides a real time accurate image of the vasculature on the surface of the skin, which helps clinicians to verify the status of the vein and reduces the risk of venipuncture.


With the help of NAVI-60 Vein Illuminator, hospital may gain the benefit of reduce costs for healthcare, improve infection control and improve patient safety and satisfaction.



Bleeding Vein      Sealed Vein      Narrow Vein      Curved Vein



NAVI-60



- 5 colors
- 4 brightness levels
- 3 projection sizes
- 2 revise modes
- 1 operation display

## + Application



### Pediatrics

Small veins "Baby fat"  
Anxiety of the children  
and their families



### Oncology

Fewer options for successful  
access because of repeated  
infusions and procedures  
and medication



### Vascular Surgery

Sclerotherapy for treating  
spider telangiectasia  
Vein stripping treatment



### Laboratory

Quick access to the vein in laboratory  
when collecting blood samples  
increases patients' satisfaction

### Portable

Hand held type Vein Illuminator, easy to carry,  
no skin contact to patient.

### Depth indicator

The depth of the vein will be detected  
and shown digitally.

### Colorful

Great image with 5 different colors optional,  
suitable for different types of skin.

