

Progreat™

Micro Catheter System

Description

[Ease of Use](#) [Description](#) [Characteristics](#) [General specifications](#) [Item specifications](#) [References](#)

Navigate distally with ease, embolize with reliability¹

Terumo offers the Progreat microcatheter system designed for therapeutic embolization and angiography in peripheral vessels. It gives control with ease of use, kink resistance, and distal flexibility in your daily practice.

Characteristics

Ease of Use

Coaxial system

Terumo's coaxial system allows fast and smooth procedures through wire integration and easy microcatheter preparation

Kink Resistance


Tungsten spiral coil & three layer structure

The tungsten spiral coil in its three layer structure protects inner lumen integrity for smooth delivery of all types of embolic materials²

Distal Flexibility

M-Coat™ & varied coil pitch

The varied pitch of the tungsten coil with Terumo's hydrophilic coating provides flexibility thus notably increasing its distal selectivity in tortuous and narrow vessels².



Progreat 2.0 Fr/2.4 Fr/2.8 Fr microcatheter (without guidewire)



Progreat 2.4 Fr/2.7 Fr/2.8 Fr microcatheter system (coaxial type with catheter and guidewire)



General specifications

Catheter inner layer	PTFE
Catheter medium layer	Radiopaque tungsten coil
Catheter outer layer	Polyester elastomer and pigment
Hydrophilic coating	M Coat polymer

MC-PP24131ZB	Progreat 2.4 Fr Coaxial	2.4 Fr / 0.80 mm	130 cm	0.022" / 0.57 mm	Straight	0.018"	0.018" / 0.46 mm	1	750 psi / 5171 kPa	0.43 ml	70 degree
MC-PC2413	Progreat 2.4 Fr Coaxial	2.4 Fr / 0.80 mm	130 cm	0.022" / 0.57 mm	Straight (shapeable)	N/A	0.018" / 0.46 mm	0	750 psi / 5171 kPa	0.43 ml	N/A
MC-PB2413	Progreat 2.4 Fr	2.4 Fr / 0.80 mm	130 cm	0.022" / 0.57 mm	Straight (shapeable)	N/A	0.018" / 0.46 mm	1	750 psi / 5171 kPa	0.43 ml	N/A
MC-PB2413ZKA	Progreat 2.4 Fr	2.4 Fr / 0.80 mm	130 cm	0.022" / 0.57 mm	Angled	N/A	0.018" / 0.46 mm	1	750 psi / 5171 kPa	0.43 ml	N/A
MC-PB2413ZKC	Progreat 2.4 Fr	2.4 Fr / 0.80 mm	130 cm	0.022" / 0.57 mm	J Curve	N/A	0.018" / 0.46 mm	1	750 psi / 5171 kPa	0.43 ml	N/A
MC-PX2413	Progreat 2.4 Fr	2.4 Fr / 0.80 mm	130 cm	0.022" / 0.57 mm	90° Angle	N/A	0.018" / 0.46 mm	0	750 psi / 5171 kPa	0.43 ml	N/A
MC-PP24151ZV	Progreat 2.4 Fr Coaxial	2.4 Fr / 0.80 mm	150 cm	0.022" / 0.57 mm	Straight	0.018"	0.018" / 0.46 mm	2	750 psi / 5171 kPa	0.47 ml	70 degree
MC-PP24151ZW	Progreat 2.4 Fr Coaxial	2.4 Fr / 0.80 mm	150 cm	0.022" / 0.57 mm	Straight	0.018"	0.018" / 0.46 mm	2	750 psi / 5171 kPa	0.47 ml	Double angled
MC-PC2415	Progreat 2.4 Fr	2.4 Fr / 0.80 mm	150cm	0.022" / 0.57 mm	Straight (shapeable)	N/A	0.018" / 0.46 mm	0	750 psi / 5171 kPa	0.47 ml	N/A
MC-PV2415Y	Progreat 2.4 Fr	2.4 Fr / 0.80 mm	150 cm	0.022" / 0.57 mm	Straight (shapeable)	N/A	0.018" / 0.46 mm	2 (3 cm proximal to distal)	750 psi / 5171 kPa	0.47 ml	N/A
MC-PV2415ZKA	Progreat 2.4 Fr	2.4 Fr / 0.80 mm	150 cm	0.022" / 0.57 mm	Angled	N/A	0.018" / 0.46 mm	2	750 psi / 5171 kPa	0.47 ml	N/A
MC-PV2415ZKC	Progreat 2.4 Fr	2.4 Fr / 0.80 mm	150 cm	0.022" / 0.57 mm	J Curve	N/A	0.018" / 0.46 mm	2	750 psi / 5171 kPa	0.47 ml	N/A
MC-PP27111	Progreat 2.7 Fr Coaxial	2.7 Fr / 0.90 mm	110 cm	0.025" / 0.65 mm	Straight (shapeable)	0.021"	0.021" / 0.53 mm	0	750 psi / 5171 kPa	0.46 ml	70 degree
MC-PE27115	Progreat 2.7 Fr Coaxial	2.7 Fr / 0.90 mm	110 cm	0.025" / 0.65 mm	Straight (shapeable)	0.021"	0.021" / 0.53 mm	0	750 psi / 5171 kPa	0.46 ml	70 degree
MC-PE27111	Progreat 2.7 Fr Coaxial	2.7 Fr / 0.90 mm	110 cm	0.025" / 0.65 mm	Straight (shapeable)	0.021"	0.021" / 0.53 mm	0	750 psi / 5171 kPa	0.46 ml	70 degree
MC-PP27131	Progreat 2.7 Fr Coaxial	2.7 Fr / 0.90 mm	130 cm	0.025" / 0.65 mm	Straight (shapeable)	0.021"	0.021" / 0.53 mm	0	750 psi / 5171 kPa	0.53 ml	70 degree
MC- PE27131	Progreat 2.7 Fr Coaxial	2.7 Fr / 0.90 mm	130 cm	0.025" / 0.65 mm	Straight (shapeable)	0.021"	0.021" / 0.53 mm	0	750 psi / 5171 kPa	0.53 ml	70 degree
MC-PC2715	Progreat 2.7 Fr	2.7 Fr / 0.90 mm	150 cm	0.025" / 0.65 mm	Straight (shapeable)	N/A	0.021" / 0.53 mm	0	750 psi / 5171 kPa	0.59 ml	N/A
MC-PP27151	Progreat 2.7 Fr Coaxial	2.7 Fr / 0.90 mm	150 cm	0.025" / 0.65 mm	Straight (shapeable)	0.021"	0.021" / 0.53 mm	0	750 psi / 5171 kPa	0.59 ml	70 degree
MC-PC2813	Progreat 2.8 Fr	2.8 Fr / 0.93 mm	130 cm	0.027" / 0.70 mm	Straight (shapeable)	N/A	0.021" / 0.53 mm	0	900 psi / 6205 kPa	0.59 ml	N/A
MC-PE28131	Progreat 2.8 Fr Coaxial	2.8 Fr / 0.93 mm	130 cm	0.027" / 0.70 mm	Straight (shapeable)	N/A	0.021" / 0.53 mm	0	900 psi / 6205 kPa	0.59 ml	N/A
MC-PP28131	Progreat 2.8 Fr Coaxial	2.8 Fr / 0.93 mm	130 cm	0.027" / 0.70 mm	Straight (shapeable)	N/A	0.021" / 0.53 mm	0	900 psi / 6205 kPa	0.59 ml	N/A
MC-PE28131ZB	Progreat 2.8 Fr Coaxial	2.8 Fr / 0.93 mm	130 cm	0.027" / 0.70 mm	Straight (shapeable)	0.021"	0.021" / 0.53 mm	1	900 psi / 6205 kPa	0.59 ml	70 degree
MC-PV2815Y	Progreat 2.8 Fr	2.8 Fr / 0.93 mm	150 cm	0.027" / 0.70 mm	Straight (shapeable)	N/A	0.021" / 0.53 mm	2 (3 cm proximal to distal)	900 psi / 6205 kPa	0.66 ml	N/A
MC-PE28151ZV	Progreat 2.8 Fr Coaxial	2.8 Fr / 0.93 mm	150 cm	0.027" / 0.70 mm	Straight (shapeable)	0.021"	0.021" / 0.53 mm	2 (3 cm proximal to distal)	900 psi / 6205 kPa	0.66 ml	N/A

 Enlarge table

Please quote above item reference codes when placing an order.

Please contact your Terumo local sales representative for more information.

References

¹ Progreat is considered as being reliable due to the kink resistance results obtained during bench testing² and Progreat is on the market for more than 15 years

² Bench testing done by Terumo Corporation. Data on file (TC-Pro107). Bench test data may not necessarily be indicative of clinical performance.

Do you need more info about this product?

Our sales representatives are here for you.

Contact us

On this page: [Description](#)[Characteristics](#)[General specifications](#)[Item specifications](#)[References](#)

Related products




Glidesheath Slender™
Introducer



Radifocus™ Guide Wire M Non-Vascular
Guide Wire



TR Band™
Radial Artery Compression Device



FemoSeal™
Vascular Closure System

Didn't find your product?

Consult our products in our [Product Index](#)