FINEÅEDIX



(1) Multi Functions

- ClearCoajet is a combination of Injection catheter and Hemostat.
- Marking, Injection, Coagulation are all possible with this device for EMR, ESD procedure.


(2) Back-cut Needle design

• It has advantages to Insertion&injection saline because the shape of Needle is sharp and I.D(Inner Diameter) is large.

Marking	Injection without exchanging device	Coagulation

SPECIFICATION

Model	Needle 0.D	Needle Thickness	Needle Length	Sheath 0.D	Working Length	Number of Use	Validity of Sterilzation
FM-CJ0002(1800) FM-CJ0002(2300)	23G	0.1	4mm	¢2.4mm	1,800mm 2,300mm	Single Use Only	3 years



SPECIFICATION

Model	Opening width	Sheath 0.D	Type of Rotation	Working Length	Number of Use	Validity of Sterilzation
FM-HG0001(1650) FM-HG0001(1820) FM-HG0001(2000)	5mm	¢2.3mm	Rotational	1,650mm 1,820mm 2,000mm	Single Use Only	3 years

Imagine What a Woderful difference You have 7