

Traditional Breast forms parameters

Silicone density - approx. 0,98 g/cm³

Product	Size	Weight, g	Volume, ml - calculated from density and weight
E-Supreme, code 509	1	86	87,7
	2	114	116,3
	3	146	148,9
	4	188	191,8
	5	250	255,1
	6	286	291,8
	7	327	333,6
	8	449	458,1
	9	572	583,6
	10	626	638,7
	11	744	759,1
	12	800	816,3
	13	999	1019,3
	14	1135	1158,1
	15	1346	1373,4
	16	1365	1392,8
	17	1385	1413,2