

Sales only to commercial customers. No sales to consumers within the meaning of §13 BGB.



Zoom in

Pressure Infusion Cuffs, reusable

The safe press can be used for liquid bags up to 500 resp. 1.000 ml filling volume. It consists of one cuff with integrated bag holder, one connecting tube with manometer and one pumpball with outlet valve.

The cuff is used for pressure-supported infusions and transfusions.

The cuff is reusable.

Allergenic potential:

Latex-free

Sterile:

No

Country availability:

Please note: not all products are available in all countries. Feel free to contact our international sales team about any questions relating to our sales and services.

Jutta Dorenberg

Director International Sales

Tel.: +49 2236 3913 - 405

eMail: jutta.dorenberg@dahlhausen.de

shorten text

Article no.	Capacity	minimal quantity
4301800013	500 ml	1
4301900013	1.000 ml	1