

## Heat sealing Sterilization Flat Roll



Apply to hospital, clinic and laboratory's sterilization; and also apply to the disinfection of beauty products or family high-temperature;

Made of transparent, polyester-polypropylene film coupled to medical grade paper which permit steam and gas penetration during sterilization.

Apply to STEAM sterilization, ETO sterilization and FORM sterilization.

The sharp part of equipment should be put contrary to the peel side to ensure safety use;

Using sealing machine to seal and mark the sterilization time, contents and usage etc.

The clear area with temperature below 25°C and humidity below 60% is recommended.

Sterility preservation period – 1 year;

The product corresponds to the standard of ISO11140-1 in conformity with EN868-5.



Cat No.	Size (mm x m)	Qty/Case (roll)
212401	Heat Sealing Sterilization Flat Roll 50x200	12
212402	Heat Sealing Sterilization Flat Roll 55x200	12
212403	Heat Sealing Sterilization Flat Roll 75x200	8
212404	Heat Sealing Sterilization Flat Roll 100x200	6
212405	Heat Sealing Sterilization Flat Roll 125x200	4
212406	Heat Sealing Sterilization Flat Roll 150x200	4
212407	Heat Sealing Sterilization Flat Roll 200x200	2
212408	Heat Sealing Sterilization Flat Roll 250x200	2
212409	Heat Sealing Sterilization Flat Roll 300x200	2
212410	Heat Sealing Sterilization Flat Roll 350x200	1
212411	Heat Sealing Sterilization Flat Roll 400x200	1
212412	Heat Sealing Sterilization Flat Roll 500x200	1