i-gel[®] supraglottic airway



For use in anaesthesia and resuscitation

i-gel® is a truly unique, single-use, latex and PVC free second generation supraglottic airway management device that represents the culmination of years of extensive research and development. Made from a medical grade thermoplastic elastomer, i-gel® has been designed to create a non-inflatable, anatomical seal of the pharyngeal, laryngeal and perilaryngeal structures whilst avoiding compression trauma.

The i-gel[®] supraglottic airway is supplied in a fully recyclable protective cradle or cage pack. This unique packaging protects the i-gel® in transit and ensures that it maintains its anatomical shape. i-gel® is available in seven sizes.



Gastric channel 15mm connector Reliable connection to any standard The i-gel® incorporates a gastric channel catheter mount or connection (except size 1). It provides an early warning of regurgitation, allows for the passing of a gastric tube to empty Proximal end of gastric channel stomach contents and facilitates venting **Clearly displayed product information** Integral bite block For quick easy reference. Includes Reduces the possibility of airway -gel confirmation of size and weight guidance channel occlusion Position guide (adult sizes only) **Buccal cavity stabiliser** Easy confirmation of optimum Aids insertion and eliminates the insertion depth potential for rotation **Epiglottic rest** Reduces the possibility of epiglottic 'down folding' and airway obstruction

The non-inflatable cuff

Made from a unique soft gel-like material allowing ease of insertion and reduced trauma

Distal end of gastric channel

Code	Description	Colour	Size		Weight	Box Qty.
8205000	i-gel [®] , supraglottic airway	•	5	Large adult	90+kg	25
8204000	i-gel®, supraglottic airway		4	Medium adult	50-90kg	25
8203000	i-gel®, supraglottic airway	•	3	Small adult	30–60kg	25
8225000	i-gel®, supraglottic airway	\bigcirc	2.5	Large paediatric	25–35kg	10
8202000	i-gel®, supraglottic airway		2	Small paediatric	10-25kg	10
8215000	i-gel®, supraglottic airway		1.5	Infant	5–12kg	10
8201000	i-gel [®] , supraglottic airway		1	Neonate	2–5kg	10

S Sterile

