


# medbar<sup>®</sup>

## Limb Holder



REF 238



CE 2292



STERILE EO

### Intended Use

Limb Holders are used for fastening the patient's hands and feet having a surgery or confined to bed who might hurt himself or environment unconsciously. Limb Holders are sterile disposable devices with single use.

### Indications

Limb Holders are used for fastening the patient's hands and feet where lying, to prevent a possible damage may be caused by patient unconsciously.

### Quality and Legision

- CE marked according to 93/42/EEC Medical Device Directive and classified as Class I sterile medical device.
- Manufactured under EN ISO 13485:2012 Quality Management System Standard.
- Notified Body: UDEM Ltd. AŞ. (2292)
- EC Certification No: M.2016.106.7000

### Biocompability

- Biological evaluation of Limb Holders done according to ISO 10993-1.
- The products do not incorporate a substance of a human blood derivative, animal orginated tissues,phthalates, medicinal product, latex, radioactive materials and electromagnetic waves.

### Sterilization

- Limb Holders are sterilized with Ethylene oxide. Validated for ISO 11135 standard.



REF NO	PRODUCT NAME	FEATURES	STERILIZATION	PACKAGING UNITS
23801	Limb Holder	• Adult Hand.	ETO	25 pcs / inner bag 250 pcs / outer box Box Dimensions (cm): 43*37*56 Box Weight: 7,00 kg
23802	Limb Holder	• Adult Foot.	ETO	25 pcs / inner bag 250 pcs / outer box Box Dimensions (cm): 43*37*56 Box Weight: 7,60 kg
23803	Limb Holder	• Kids.	ETO	100 pcs / inner bag 500 pcs / outer box Box Dimensions (cm): 39*43*41 Box Weight: 6,50 kg
23804	Limb Holder	• Baby.	ETO	100 pcs / inner bag 500 pcs / outer box Box Dimensions (cm): 39*43*41 Box Weight: 6,50 kg

