

## Weight & Size overview

### Front aprons

| Article code | Sizes       | Weight 0.35 mm Leq | Weight 0.5 mm Leq | Chest size | Length |
|--------------|-------------|--------------------|-------------------|------------|--------|
| Front-S      | Small       | 2.7 kg             | 3.6 kg            | Up to 90   | 100    |
| Front-M      | Medium      | 3.1 kg             | 4.1 kg            | Up to 100  | 100    |
| Front-L      | Large       | 3.4 kg             | 4.5 kg            | Up to 110  | 100    |
| Front-XL     | Extra-large | 3.7 kg             | 4.9 kg            | Up to 120  | 105    |

### UFront aprons

| Article code | Sizes       | Weight 0.35 mm Leq | Weight 0.5 mm Leq | Chest size | Length |
|--------------|-------------|--------------------|-------------------|------------|--------|
| UFront-S     | Small       | 2.9 kg             | 3.9 kg            | Up to 90   | 110*   |
| UFront-M     | Medium      | 3.3 kg             | 4.4 kg            | Up to 100  | 115*   |
| UFront-L     | Large       | 3.6 kg             | 4.8 kg            | Up to 110  | 120*   |
| UFront-XL    | Extra-large | 3.9 kg             | 5.2 kg            | Up to 120  | 120*   |

\* Longer than our standard front apron in able to work with it while sitting down.

### Coat aprons

| Article code | Sizes male   | Weight 0.5/0.25 mm Leq | Weight 0.7/0.35 mm Leq | Chest size | Hip     | Length |
|--------------|--------------|------------------------|------------------------|------------|---------|--------|
| Coat-A       | Small        | 4.0 kg                 | 5.6 kg                 | Up to 90   | 95-100  | 100    |
| Coat-B       | Medium       | 4.3 kg                 | 6.0 kg                 | Up to 110  | 110-115 | 100    |
| Coat-C       | Large        | 4.4 kg                 | 6.2 g                  | Up to 115  | 115-120 | 100    |
| Coat-CXL     | Extra-large  | 5.0 kg                 | 7.0 kg                 | Up to 120  | 120-125 | 105    |
| Article code | Sizes female | Weight 0.5/0.25 mm Leq | Weight 0.7/0.35 mm Leq | Chest size | Hip     | Length |
| Coat-D       | Small        | 3.8 kg                 | 5.3 kg                 | Up to 85   | 90-95   | 100    |
| Coat-E       | Medium       | 4.1 kg                 | 5.7 kg                 | Up to 95   | 95-100  | 100    |
| Coat-F       | Large        | 4.2 kg                 | 5.9 kg                 | Up to 105  | 100-110 | 100    |
| Coat-FXL     | Extra-large  | 4.8 kg                 | 6.7 kg                 | Up to 115  | 115-120 | 105    |

### Vest & Kilt aprons

| Article code | Sizes male   | Weight 0.5/0.25 mm Leq | Weight 0.7/0.35 mm Leq | Chest size | Hip     | Vest length | Kilt length |
|--------------|--------------|------------------------|------------------------|------------|---------|-------------|-------------|
| VK-A         | Small        | 4.7 kg                 | 6.5 kg                 | Up to 90   | 95-100  | 57          | 59          |
| VK-B         | Medium       | 6.3 kg                 | 8.8 kg                 | Up to 110  | 110-115 | 62          | 64          |
| VK-C         | Large        | 6.8 kg                 | 9.5 kg                 | Up to 115  | 115-120 | 65          | 67          |
| VK-CXL       | Extra-large  | 7.4 kg                 | 10.3 kg                | Up to 120  | 120-125 | 68          | 70          |
| Article code | Sizes female | Weight 0.5/0.25 mm Leq | Weight 0.7/0.35 mm Leq | Chest size | Hip     | Vest length | Kilt length |
| VK-D         | Small        | 4.5 kg                 | 6.3 kg                 | Up to 85   | 90-95   | 51          | 55          |
| VK-E         | Medium       | 4.8 kg                 | 6.7 kg                 | Up to 95   | 95-100  | 55          | 58          |
| VK-F         | Large        | 5.4 kg                 | 7.5 kg                 | Up to 105  | 100-110 | 56          | 61          |
| VK-FXL       | Extra-large  | 6.2 kg                 | 8.6 kg                 | Up to 115  | 115-120 | 57          | 63          |

0.5mm Leq = 6kg/m2

Weight according to EN 61331-3:2014, this list can only be used as a reference, all listed weights have a tolerance of 7%. W/S = 0.25/0.25mm Leq = 3kg/m2 / W/S = 0.35/0.35mm Leq = 4,2kg/m2



### Thyroid shields

3 different models, adjustable size using magnet closure

0.5 mm Leq

| Article code | Weight*   |
|--------------|-----------|
| VIPE         | 250 grams |
| TM-3         | 340 grams |
| TS-4         | 320 grams |

### Gonad Protection

0.5 mm Leq

| Article code | Weight*   |
|--------------|-----------|
| GP-S         | 600 grams |
| GP-M         | 1.0 kg    |
| GP-L         | 1.4 kg    |
| GP-Belt      | 750 grams |

### Head Protection

0.5 mm Leq

| Article code | Weight*   |
|--------------|-----------|
| HP-S         | 250 grams |
| HP-M         | 300 grams |
| HP-L/LX      | 350 grams |
| HP-XXL       | 400 grams |

0.5mm Leq = 6kg/m<sup>2</sup>

Weight according to EN 61331-3:2014, this list can only be used as a reference, all listed weights have a tolerance of 7%. W/S = 0.25/0.25mm Leq = 3kg/m<sup>2</sup> / W/S = 0.35/0.35mm Leq = 4,2kg/m<sup>2</sup>