


20.10.12

Dispomed Witt oHG

Am Spielacker 10 -12, D-63571 Gelnhausen
Tel.: (+49) 0 60 51 - 96 23 -0
Fax.: (+49) 0 60 51 - 96 23 -0
eMail: info@dispomed.de




Dispomed®

PERFUJECT® Pump syringes

3-part syringes 50 ml, adapted for Braun Perfusor and compatible devices

- According to DIN EN ISO 788
- Luer Lock, centric
- Pressure tested to 6 bar
- Rub firm graduation
- Cylinder and plunger made of Polypropylene
- latexfree
- pistonstop
- exact reading of the syringe- volume by clearly 1 ml graduation
- for application of photosensitive pharmaceuticals with orange(amber) or black colouring.
- Sterile
- R (Gamma radiation)




			Art.-No.	pack. unit
50 ml, without needle, transparent			22051	50 / 300
<u>50 ml, transparent, with attached needle</u>	<u>2,0 x 30mm</u>	<u>14G x 1^{1/4}"</u>	<u>22052</u>	50 / 300
50 ml, black, with attached needle	2,0 x 30mm	14G x 1 ^{1/4} "	22053	50 / 300
50 ml, orange, with attached needle	1,8 x 25mm	15G x 1"	22063	50 / 300

Black coloured syringes have a red graduation on the plunger

Dispomed

Lotus 12

Dispomed Witt oHG
Am Spielacker 10 -12, D-63571 Gelnhausen
Tel.: (+49) 0 60 51 - 96 23 -0
Fax.: (+49) 0 60 51 - 96 23 -0
eMail: info@dispomed.de




Dispomed®

ECOJECT® PLUS 3-part syringe

ECOJECT® PLUS

3-part syringe LUER LOCK with elastomere piston

- accurate guidance of the plunger
- precise dosing
- piston plug with twice sealing ring
- rub-firm graduation
- Latex and PVC-free
- sterile EO (Ethylene oxide)



ECOJECT® PLUS Luer Lock

	Art.-No.	pack. unit
ECOJECT® PLUS 2,5 ml	21004	100 / 3000
ECOJECT® PLUS 5 ml	21005	100 / 2000
ECOJECT® PLUS 10 ml	21010	100 / 1200
ECOJECT® PLUS 20 ml	<u>21020</u>	50 / 1000

Asier