



8403DXC

C-MAC® Video Laryngoscope MIL #0,

CMOS technology, with MILLER laryngoscope blade, size 0, with C-MAC® system interface, documentation of images and video sequences via the BlueButton, for use with C-MAC® Connecting Cable 8403X for C-MAC® Monitor 8403ZXK or C-HUB® II 20290301 as well as C-MAC® PM 8403XD



8403GXC

C-MAC® Video Laryngoscope MIL #1,

CMOS technology, with MILLER laryngoscope blade, size 1, with C-MAC® system interface, documentation of images and video sequences via the BlueButton, for use with C-MAC® Connecting Cable 8403X for C-MAC® Monitor 8403ZXK or C-HUB® II 20290301 as well as C-MAC® PM 8403XD



8403EXC

C-MAC® Video Laryngoscope MAC #0,

CMOS technology, with MACINTOSH laryngoscope blade, size 0, with C-MAC® system interface, documentation of images and video sequences via the BlueButton, for use with C-MAC® Connecting Cable 8403X for C-MAC® Monitor 8403ZXK or C-HUB® II 20290301 as well as C-MAC® PM 8403XD



8403KXC

C-MAC® Video Laryngoscope MAC #2,

CMOS technology, with MACINTOSH laryngoscope blade, size 2, with C-MAC® system interface, documentation of images and video sequences via the BlueButton, for use with C-MAC® Connecting Cable 8403X for C-MAC® Monitor 8403ZXK or C-HUB® II 20290301 as well as C-MAC® PM 8403XD



8403HXP

C-MAC® Pediatric Video Laryngoscope D-BLADE,

CMOS technology, for difficult intubation in children, with C-MAC® system interface, documentation of images and video sequences via the BlueButton, for use with C-MAC® Connecting Cable 8403X for C-MAC® Monitor 8403ZXK or C-HUB® II 20290301 as well as C-MAC® PM 8403XD