

3.1. Table controls

For table controls see section [Patient table control](#) section.

4. Multimodality patient table

If you want to purchase a CT or PET-CT to the SPECT gantry later, then practical to order carbon-fiber table to the system.

The carbon fiber dual section patient table can be extended extremely long, ensuring safe support for obese patients as well.

- Intelligent CAN connected industry standard interface to gantry.
- Dual tabletop.
- Motorized vertical movements.
- Motorized whole body motion.
- Low attenuation (< 9% at 140 keV) carbon fiber pallet.
- Horizontal moving range is 388.5 cm (max. scan length for WB and WB SPECT is 245 cm with extender, for SPECT-CT and Helical CT is 205 cm, Localizer length is 202 cm).
- Height of the patient pallet is variable between 46-83 cm, optionally 95 cm.
- 229 kg (505 lbs) patient weight.