

## Configuratie DERMATOME Skin Graft Mesher BES Battery

1. K6-252 CONTAINER 3/4 STANDARD MODEL - GREY 1 BOTTOM NON-PERFORATED - LID PERFORATED
2. K2-740 SINGLE USE PAPER FILTER WITH INDICATOR 1 MARKED - (PACK OF 100 PCS)
3. K2-660 SECURITY SEAL - RED 1 (PACK OF 100 PCS) - WITHOUT INDICATOR
4. K2-672 PAPER LABEL 1 WITH INDICATOR - (PACK OF 100 PCS)
5. K2-673 PAPER LABEL 1 WITH INDICATOR - (PACK OF 100 PCS)
6. M11-500 SKIN GRAFT MESHER
7. M11-600 SKIN GRAFT CARRIER 1,5:1 (BOX OF 20)
8. M11-605 SKIN GRAFT CARRIER 3:1 (BOX OF 20)