

Occlusion Balloon Catheters

The New Standard of Excellence, Redefining Deliverability, Versatility, and Control

Scepter	Scepter C Compliant Balloon Catheter				
Product Code	Description	Balloon Diameter (mm)	Balloon Length (mm)	Distal Tip Length (mm)	
1 balloon catheter, 1 introducer sheath, 1 shaping mandrel per box					
BC0410C	Compliant Occlusion Balloon 4.0mm x 10mm	4.0	10	5	
BC0415C	Compliant Occlusion Balloon 4.0mm x 15mm	4.0	15	5	
BC0420C	Compliant Occlusion Balloon 4.0mm x 20mm	4.0	20	5	
ScepterXC X-tra Compliant Balloon Catheter					
Product Code	Description	Balloon Diameter (mm)	Balloon Length (mm)	Distal Tip Length (mm)	
1 balloon catheter, 1 introducer sheath, 1 shaping mandrel per box					
BC0411XC	X-tra Compliant Occlusion Balloon 4.0mm x 11mm	4.0	11	5	

reatures	Denetits
14 Wire Compatible	Provides more choices of guidewire (Traxcess® docking wire compatible) Improves trackability and stability with better performance of the wire
Hydrophilic Coating on the Balloon	Significantly improves trackability Enables use of longer balloons without compromising trackability
Soft & Long Tip	Improves balloon stability and trackability
Low Profile	Ensures compatibility with double catheter technique in 6F guide Improves trackability
Distal Tip Marker	Positive identification of distal tip location

Ronofite

MicroVention, Inc. Worldwide Headquarters 1311 Valencia Avenue Tustin, CA 92780 USA MicroVention UK Limited MicroVention Europe, S.A.R.L. MicroVention Deutschland GmbH Web



PH +44 (0) 191 258 6777 PH +33 (1) 39 21 77 46 PH +49 211 210 798-0 **microvention.com**



Easturas

MICROVENTION and Traxcess are registered trademarks of MicroVention, Inc. Scepter C and Scepter XC are trademarks of MicroVention, Inc. • Scientific and clinical data related to this document are on file at MicroVention, Inc. Refer to Instructions for Use for additional information. ©2013 MicroVention, Inc. MM228(I) Rev.C 1/13