

Disposable Endoscopic linear cutting stapler and accessories



Indication

1. Suitable for open or endoscopic surgery
2. Resection, transection and stapling of lung, bronchial tissue, stomach and intestine.

Characters

- ❖ Strengthened anvil, better grip power; slim shape, stronger passability.
- ❖ The stepped staple cartridge is gradually clamped, making it easier to press the tissue to the desired height.
- ❖ High and low staples, adapt to the natural state of the organization; tight inside and loose outside, reducing the occurrence of anastomotic leakage.
- ❖ The drop-shaped slot design makes it easier to molding.
- ❖ 11 angle adjustments, maximum 45° on each side.
- ❖ The new enlarged unloading button design makes the operation more convenient.
- ❖ Ring handle, one-handed operation; control jaws to clamp and release tissue.

Stapler No.	Cartridge No.	Cartridge Color	Total Length	Suturing Length	Suturing Thickness			Accessory type	
					First Row	Second Row	Third Row		
	Lot 2.2	JJIIID301	Brown	202	30	2.0	2.5	3.0	Rotating type
		JJIIID302	purple	202	30	3.0	3.5	4.0	Rotating type
		JJIIID303	black	202	30	4.0	4.5	5.0	Rotating type
JQIIA(S)	Lot 2.3	JJIIID451	Brown	217	45	2.0	2.5	3.0	Rotating type
JQIIB(Short)		JJIIID452	purple	217	45	3.0	3.5	4.0	Rotating type
JQIIC(Long)		JJIIID453	black	217	45	4.0	4.5	5.0	Rotating type
	Lot 2.4	JJIIID601	Brown	232	60	2.0	2.5	3.0	Rotating type
	Lot 2.7	JJIIID602	purple	232	60	3.0	3.5	4.0	Rotating type
		JJIIID603	black	232	60	4.0	4.5	5.0	Rotating type

