# **ActionBLU**:

(Trypan Blue Ophthalmic Solution 0.06% w/v)



## **Description:**

ActionBlu 0.06% is a dye which safely stains the anterior lens capsule & is mainly used in cataract surgery. Special feature of Optiblu is that is stains the anterior capsule without affecting the corneal endothelium.

# Indication:

ActionBlu is indicated for facilitating surgical procedure for cataract extraction (IOL Surgery and phacosurgery) and as a staining agent.

## **Contraindication:**

ActionBlu is contraindicated in patients with hypersensitivity to the dye or any of its components.

# Warning & Precautions:

- Do not use the product after the Expiry date
- Do not use on children
- For single use only
- Not for injection

#### Storage:

Store at a temperature between 15° Celsius and 35° Celsius and relative humidity less than 65% protect from freezing.

Do not expose containers to direct sunlight

# Supply:

ActionBlu is a sterile solution is available in 1ml with pouch in pack of Ten and ActionBlu PFS 1ml in sterile glass syringe.

Another Quality Product from:

Action

# Action Medical Pvt Ltd.

(An ISO 9001:2008 & ISO 13485:2003 Certified Company) Gat No 528, Koregaon-Bhima, Near SBI, Pune-Nagar Road Taluka Shirur, Dist Pune Pin 412216, INDIA Tel: 0091-2137-253617 / 18. Fax: 0091-2137-253620. E-mail : <u>action@actionmedihealth.com</u> Website: www.actionmedihealth.com