

Infusion/Transfusion

Pressure Infusion Cuffs



Pressure Infusion Cuffs

The safe press can be used for liquid bags up to 500 resp. 1.000 ml filling volume. It consists of one cuff with integrated bag holder, one connecting tube with manometer and one pumpball. The cuff is used for pressure-supported infusions and transfusions.

Sterile: No

REF	Capacity	PU
43.018.00.012	500 ml	20 pcs
43.019.00.012	1.000 ml	20 pcs



Pressure Infusion Cuffs, reusable

The safe press can be used for liquid bags up to 500 resp. 1.000 ml filling volume. It consists of one cuff with integrated bag holder, one connecting tube with manometer and one pumpball with outlet valve. The cuff is used for pressure-supported infusions and transfusions. The cuff is reusable.

Allergenic potential: Latex-free

REF	Capacity	PU
43.018.00.013	500 ml	6 pcs
43.019.00.013	1.000 ml	6 pcs