










# PRE-ASSEMBLED EEG HEAD CAPS WITH STANDARD-THICKNESS ELECTRODES AND WIRES LOCATED EXTERNALLY **REUSABLE**

**Packaging: box containing 1 piece.**

CODE	WIRES	SIZE	Number of electrodes	DESCRIPTION	COLOUR
<b>CUEE60S</b>	external	58-62 cm	21, standard thickness	Extra large, adult	Blue 
<b>CUEE56S</b>	external	54-58 cm	21, standard thickness	Large, adult	Red 
<b>CUEE52S</b>	external	50-54 cm	21, standard thickness	Medium-small, adult	Yellow 
<b>CUEE48S</b>	external	46-50 cm	21, standard thickness	Large, child	Green 
<b>CUEE44S</b>	external	42-46 cm	21, standard thickness	Medium, child	L. Blue 
<b>CUEE41S</b>	external	38-42 cm	21, standard thickness	Small, child	Pink 
<b>CUEE38S</b>	external	34-38 cm	11, standard thickness	Large, neonatal	Brown 
<b>CUEE34S</b>	external	32-34 cm	11, standard thickness	Medium, neonatal	Orange 
<b>CUEE32S</b>	external	30-32 cm	11, standard thickness	Small, neonatal	White 



**CUSTOMIZED CONNECTION CABLES ARE AVAILABLE UPON REQUEST.**

YOU CAN ORDER EEG HEAD CAPS WITH DIRECT CONNECTION TO THE EEG BOX.  
JUST REPLACE THE PRODUCT CODE'S LAST LETTER (**S**) WITH THE LETTER STANDING FOR THE  
DESIRED CONNECTION. SEE THE CHART BELOW.

<b>A</b>	INTERNATIONAL		<b>N</b>	NEUROSOFT
<b>M</b>	MICROMED		<b>XX</b>	CUSTOM-MADE

# EAR ELECTRODE

## REUSABLE

This electrode with elastic clamp must be placed on the ear lobe.  
The electrode is made up of a plastic body and an Ag/AgCl recording area.

**Packaging: box containing 1 piece.**

CODE	DIAMETER	CABLE LENGTH	INSTRUMENT END
<b><i>EARPLD1-10</i></b>	10 mm	1.0 m	Ø 1.5 mm touchproof DIN42802

