




GEISTER®

Neuro and Spine Surgery
Orthopaedic Surgery





Lumbar Retractor System

quick connect


Fig.	Art.-No.	Description
A	29-3605	Transversal Lumbar Retractor, hinged, anterior titanium hinges



Introducer handle for blades

Can also be used as a hand-held retractor

Fig.	Art.-No.	Description
B	29-3615.FL	black PPSU handle, 195 mm / 7 1/2" Luer-Lock and locking mechanism, with Luer-Lock seal cap, silicone



Depth Gauge, PEEK nature


to determine the appropriate blade length

Fig.	Art.-No.	Description
C	29-3609	160 mm, measuring range 100 mm



Lumbar Retractor System – SpineControl™ (quick connect)


GEISTER®



Lumbar Retractor Set

complete in tray with lid
quick connect, Titanium

29-3604 Consisting of			
Art.-No.	Qty.	Art.-No.	Qty.
29-3605	1	29-3651.65	2
29-3615.FL	2	29-3651.75	2
29-3609	1	29-3652.35	2
29-3650.35	2	29-3652.40	2
29-3650.40	2	29-3652.45	2
29-3650.45	2	29-3652.55	2
29-3650.55	2	29-3652.65	2
29-3650.65	2	29-3652.75	2
29-3651.35	2	29-3652.85	2
29-3651.40	2	29-3655.01	1
29-3651.45	2	29-3655.02	1
29-3651.55	2	77-7286	1



Lumbar Blades

20 mm width, Titanium, quick-connect

Fig.	Art.-No.	Length mm
D	29-3650.35	35
D	29-3650.40	40
D	29-3650.45	45
D	29-3650.55	55
D	29-3650.65	65

40 mm width, Titanium, quick-connect

Fig.	Art.-No.	Length mm
E	29-3651.35	35
E	29-3651.40	40
E	29-3651.45	45
E	29-3651.55	55
E	29-3651.65	65
E	29-3651.75	75

60 mm width, Titanium, quick-connect

Fig.	Art.-No.	Length mm
F	29-3652.35	35
F	29-3652.40	40
F	29-3652.45	45
F	29-3652.55	55
F	29-3652.65	65
F	29-3652.75	75
F	29-3652.85	85

Tray

for SpineControl™ Lumbar System

Art.-No.	Description
77-7286	1/1 Type, with lid, 485×253×70 mm, perforated stainless steel



Hemi-Laminectomy Hook

Titanium, quick-connect, blue

Fig.	Art.-No.	Description
G	29-3655.01	small
G	29-3655.02	large