



pelicangio

Syringes for contrast media injectors

- Wide range of syringes adaptable for Medtronic, Medrad and Covidien injectors



Features & benefits

RESISTANT TO 300 PSI (21 BARS) OR 1200 PSI (83 BARS)

- For low or high pressure injections

CAPACITY

- From 130ml to 200ml

References

Medrad injectors

References	Volume (ml)	Application	Injectors	Accessories	Packaging
0137NB	130	Angio.	Mark III & Mark IV	Fill tube	50/box
0137NN	150	Angio.	MarkV,V+ & Provis	'Short nose' version + fixed Luer-lock + fill tube	
0137ND	200	Angio.	MarkV,V+ & Provis	Fill tube	
0137NH	200	Angio.	MarkV,V+ & Provis	-	
0137RA	200	CT	MarkV,V+ & Provis	Fill tube + coiled 150cm connection line	
0137RB	200	CT	MarkV,V+ & Provis	Coiled 150cm connection line + spike	
0137RC	130	Angio.	Mark III & IV	Coiled 150cm connection line + spike	

Covidien injectors

References	Volume (ml)	Application	Injectors	Accessories	Packaging
0137NB	130	CT	Optistar	Fill tube	50/box
0137NE	150	CT & angio.	Angiomat 6000 & Illumena	Fill tube	
0137RC	130	Angio	Optistar	coiled 150cm connection line + spike	

Medtron injectors

References	Volume (ml)	Application	Injectors	Accessories	Packaging
0137NB	130	CT, MRT & angio.	Injektron 82CT, CT2, MRT & 82HP	Fill tube	50/box
0137ND	200	CT, MRT & angio.	Injektron 82CT, CT2, MRT & 82HP	Fill tube	
0137NS	Kit 2 x 200	CT	Injektron 82CT & CT2	2 fill tubes + 1 x 150 cm straight line	20/box
0137RA	200	CT	Injektron 82CT & CT2	Fill tube + coiled 150cm connection line	50/box
0137RB	200	CT	Injektron 82CT & CT2	Coiled 150cm connection line + spike	
0137RC	130	Angio.	Injektron 82HP	Coiled 150cm connection line + spike	
0137NH	200	Angio.	Injektron 82HP	—	50/box

Please read carefully the Instructions For Use. The CE conformity assessment has been conducted by LNE-GMED (0459).

Manufacturer:  PEROUSE MEDICAL - Route du Manoir - 60173 Ivry le Temple - France - Tel: +33 (0) 3 44 08 17 00